


TECHNICAL STUDY



Nucleation



**Cycle time
reduction**



**Reduced
costs**



Demolding




THE EFFECT OF CEVO® 2443
ON THE NUCLEATION BEHAVIOUR OF PET

Introduction

CEVO® 2443 is a proprietary formulation based on saponified, polar montan wax with a particularly high proportion of calcium montante.

In the present study the crystallisation temperature of PET with 0,5 % CEVO® 2443 has been DSC-analysed

in comparison to the non additivated polymer. The aim of the study was to demonstrate a nucleation effect of CEVO® 2443 by increasing the PET crystallization temperature.



	Target value	Unit	Method
Physical form	yellowish compactate		
Acid number	< 10	[mg KOH/g]	ISO 2114
Melting range	ca. 120	[°C]	DSC
Na content	min. 5	[%]	RFA

Table 1: Characteristics of CEVO® 2443

Test material

1. PET: Lighter C93, unfilled (Equipolymers)

2. CEVO® 2443 (VOELPKER)

Pretreatment of PET: dried at 120 °C / 2,5 h

Test equipment

Injection molding: Arburg Allrounder 320 C 500-170

Temperature profile: 250/265/275/285 °C

Differential Scanning Calometry (DSC):
Mettler Modul: 822e/1344/500

Method: PET ASTM D 3410 40-300:

Dt 0,5 s

[1] 40 °C, 5 min

[2] 40,0..300 °C, 10,00 K/min

[3] 300°C, 5 min

[4] 300,0 °C..40 °C, -2,00 K/min

[5] 40,0 °C, 5 min

[6] 40,0 °C..300,0 °C, 10 K/min

Results




Table 2:
Shift of the DSC Crystallisation temperature of PET
→ and PET with CEVO® 2443




Table 4: Difference of the DSC Crystallisation temperature PET and PET with CEVO® 2443

Nr.	PET/PET Formulation	Solidification temperature [°C]	Solidification temperature Difference [Δ °C]
1	100 % PET	212,8	-
2	99,5 % PET + 0,5 % CEVO® 2443	228,7	+ 15,9

Table 3: Formulation and solidification temperature difference

Summary

CEVO® 2443 acts as a suitable nucleating agent in PET. It can enhance the crystallization rate of PET and provide reduced cycle times in injection molding.

CEVO® 2443 is effective at concentrations of 0.1 – 0.5 % in polyester. It exhibits additional desirable

properties for PET injection molding, for example: low volatility and favourable lubrication, mold-release and surface improving effects. CEVO® 2443 thus contributes to improved PET processing and final product performance.

Völpker Spezialprodukte GmbH
Fabrikstraße 1 | 39393 Völpke | Germany
Tel. +49 (0) 39402 962-0
Fax +49 (0) 39402 215
application@voelpker.com
www.voelpker.com

Design: www.artfaktor.de



This flyer is printed on
certified recycled paper.

ecovadis



1st Edition | 02/2025

Disclaimer

The information contained herein are believed to be accurate and reliable as of the date issued. However, we do not warrant or guarantee the accuracy or reliability, in particular not for any specific intended use by the customer. It is the responsibility of those to whom we supply our products directly or indirectly to ensure that their use of the product complies with existing regulations, laws, legislations and proprietary rights. The information given by Völpker Spezialprodukte does not exempt the customer from carrying out inspections and analyses on goods purchased.